

CONTENTS

Introduction How to Create Your Shoe - Nimco Mac Toddlers Footwear Childrens Footwear Boa Closure System Waterproof Boots Adaptions Bespoke to Measures or Cast General Terms Measurements and Lasts	de4You	3 4 6 14 50 52 54 56 57 58

New For 2020

We are proud to present the 2020/2021 Nimco Made4You Childrens Footwear Range.

What's New:

- 17 New Styles
- New colours for existing styles
 - Waterproof Boots
- New Featherweight sole unit
 - New Extra Narrow Last



Introduction

The Nimco Made4You range is now available in 39 different styles which have been subdivided into 2 different categories: Toddlers styles (covering sizes 16 - 27) & Childrens styles (covering sizes 22 - 45).

As always, most styles can be ordered in lace or velcro and many with a choice of 2 back heights: standard or continental (25mm higher)

All Nimco shoes are personally made for you, this allows us to provide you with an unmatched level of reliability. All standard complete orders will be dispatched within 12 to 14 working days.

We are also able to offer a full range of adaptions to be added to your order. The majority of adaptions that we provide are applied at the time of manufacture which increases the overall quality whilst maintaining our reliable lead time.



Create your footwear in 4 steps!

1. Style

When creating your new footwear, the first step is to select your desired style. Once you have decided, you can then select the correct size, closure and backheight.



2. Last and Stiffener

Once you have selected your style and options, you can select the required last and stiffener. These selections will determine the type of support provided by your new footwear and ensure that they fit perfectly.

Stiffener Options

Normal Stiffener



This stiffener is ideal for standard lasts and with insoles and supports.

Stability Stiffener



The solid butterfly shaped heelcap in conjunction with the high upper of a boot offer optimum support

Orthosis Stiffener



This sturdy heel cap compliments orthoses lasts as these are shaped to accommodate an AFO

Anti-Varus Heel & Toe Cap



The AV stiffeners strengthen the shoe upper in the areas required to treat Metatarsus Adductus

Last Options



All styles are available in 14 different last options: 6 Standard Lasts 3 Orthosis Lasts 5 Anti Varus Lasts

Please consult page 58 for sizes/widths



3. Outsole

Your choice of outsole will be influenced on whether a stitchdown or cemented sole best meets your needs.

PU Sole (Flexible and thin PU)








Available in Cemented and Stitchdown construction

Sport Sole (EVA midsole outsole made of Rubber)











White

Available in Cemented and Stitchdown construction

Rubber Sole (Long lasting strong rubber)





Black





Stitchdown Sole

CEMENTED Sole

Available in Cemented and Stitchdown construction


Welt Sole (Very light weight and flexible)





Rubber









Black Spot White TR Ultralight TR Ultralight

TR Sole* (Ultralight weight upper layer and TR bottom material)







*Please be aware when ordering footwear with this sole, flares and sockets won't be possible, please select one of the other sole units available.

4. Colours & Materials

You can specify which colour/material you require on any aspect of the upper. This is also an ideal solution if you would like to make a particular style suitable for use at school by selecting an appropriate colour such as black, brown or blue. Please state clearly which colours you require on your order.

Please see our website for all the colours available.

Please note: that due to the high level of personalisation, the addition of custom colours to your order will mean that they are treated as non returnable.



Toddlers







R1995/K1906 Waxed Leather Brown Combi



NM1902/R552 Waxed Leather Black Combi



NM1997/L1676 Nappa Blue Combi

Closure Lace or Velcro

Sizes 16 - 27 **H**eiqhт UK Only

Standard & Stability 6 widths

Lasts & Stiffeners **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles are



TR Sole





PU Sole















Rubber Sole













NM1902/P1977 Nappa Green Combi



NM1997/P1697 Blue Nappa Combi



NM1902/P402 Waxed Leather Black Combi

Closure Lace or Velcro

Sizes 16 - 27 **НеіqҺ**т UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

5 widths

Soles =

The following 2 soles



TR Sole





PU Sole







Sports Sole























NM1997/Y1997 Fantasy Blue Combi



NM1902/Y1915 Fantasy Purple Leather Combi



NM1902/Y1972 Fantasy Black Combi



N325/Y1925 Fantasy Pink Combi

Closure Lace or Velcro

Sizes 16 - 27 **H**eiqhт UK Only

Lasts & Stiffeners
Standard & Stability Orthosis 6 widths

3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole







Rubber Welt Sole









Y1955 Fantasy Beige





S602 Patent Black



R577 Waxed Leather Fuchsia

Options

Closure Velcro Only **Sizes** 16 - 27

НеіqҺт UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles



TR Sole





PU Sole







Sports Sole



Black Brown White

























P1672 Split Suede Blue



Y1508 Fantasy Black

Closure Lace or Velcro **Sizes** 16 - 27

Heiqhт UK Only Standard & Stability
6 widths

LASTS & STIFFENERS Orthosis 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles are



TR Sole





PU Sole





























L1601/P496 Suede Fuchsia Combi



L1601/P497 Suede Blue Combi



NM1902/Y1972 Fantasy Black Combi

Closure Lace or Velcro

Sizes 16 - 27

Неіднт UK Only Standard & Stability 6 widths

Lasts & Stiffeners **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole





Rubber Welt Sole











NM1997/TG1903 Fabric Fantasy Blue Combi



NM1902/Y1932 Fantasy Black Combi



NM1902/TG1902 Fabric Fantasy Black Combi

Closure Lace or Velcro

Sizes 16 - 27 **НеіqҺ**т UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

5 widths

Soles =

The following 2 soles



TR Sole





PU Sole













Black Brown White













Childrens









J846/M1402 Polished Leather Black Combi



R502/K1902 Nappa Black Combi



V1430/Y693 Nappa Blue Combi

Closure Lace or Velcro or BOA

Sizes 22 - 45 **H**e**iqh**т UK Only

Lasts & Stiffeners Standard & Stability Orthosis 6 widths

3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole

















L1602 Nappa Black Combi

Closure Lace or Velcro

Sizes 22 - 45 **НеіqҺ**т UK Only

Standard & Stability 6 widths

Lasts & Stiffeners **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles



TR Sole





PU Sole







Sports Sole



Black Brown White







Rubber Sole















S602 Patent Black

Closure Velcro Only

Sizes 22 - 45

Heiqhт UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles



TR Sole





PU Sole























Touchdown (Sports sole as standard)









V1420 Nappa Navy

*Please be aware that both Brown and Blue images are shown in Continental height.

Options

Closure Lace or Velcro

Sizes 22 - 45 **Неід**Інт UK & Continental

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis**

3 widths

Anti-Varus 5 widths

Soles

The following 2 soles are



TR Sole



Spot White

PU Sole

















Black Brown White













R552/Y692 Waxed Leather Black Combi



R1955/P407 Waxed Leather Brown Combi

Closure Lace or Velcro

Sizes 22 - 45 **H**eiqhт UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles



TR Sole





PU Sole















Rubber Sole



















L1672/Y1508 Fantasy Black Combi



N302/S1997 Fantasy Blue Combi

Closure Lace or Velcro

Sizes 22 - 45 Height UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles

The following 2 soles



TR Sole





PU Sole



























ASHTON (PU sole as standard)





R552/TS1976 Waxed Leather Black Combi



J846/L1602 Polished Leather Black



M1654/V1430 Full Grain Aniline Brown Combi

Options

Closure Lace or Velcro

Sizes 22 - 45 Height UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles

The following 2 soles



TR Sole





PU Sole













Black Brown White













Y1972/P402 Fantasy Black





S602/P402 Patent Black



V1430/P497 Nappa Blue



R502/N302 Nappa Black



N328/P497 Nubuck Navy Combi



S1505/P496 Patent Fuchsia With Dots *this leather will incur a surcharge

Closure Lace or Velcro or BOA

Sizes 22 - 45 **Неі**q**н**т UK & Continental

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles

The following 2 soles are



TR Sole





PU Sole

































R555/P1668 Waxed Leather Brown Combi



R502/AR1977 Nappa Black Combi



R552/P452 Waxed Leather Black Combi

Closure Lace or Velcro

Sizes 22 - 45 Heiqhт UK & Continental

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles

The following 2 soles



TR Sole





PU Sole

















Black Brown White















E20122/Y1915 Fantasy Purple Combi



P407/Y1955 Fantasy Brown Combi



N302/Y1932 Fantasy Black Combi

Closure Lace or Velcro

Sizes 22 - 45 **H**е**i**q**h**т UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

5 widths

Soles =

The following 2 soles



TR Sole





PU Sole





Sports Sole









Rubber Sole













L1601/Y1995 Fantasy Beige Combi



L1601/S1825 Fantasy Pink Combi



L1601/S602 Fantasy Black Combi



L1601/Y1632 Fantasy Fuchsia Combi

Closure Lace or Velcro or BOA

Sizes 22 - 45

Неіднт UK Only

Lasts & Stiffeners Standard & Stability 6 widths

Orthosis 3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole















FANTASY (Sports sole as standard)





Y1972/Y1634 Fantasy Black Combi



N325/Y1663 Nubuck Pink Combi



Y1925/N325 Fantasy Pink Combi

Options

Closure Velcro Only

Sizes 22 - 45 **НеіqІнт** UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles are



TR Sole





PU Sole













Black Brown White











JASMIN (PU sole as standard)



E20122/Y1632 Nappa Purple Combi



Y1932/L1602 Fantasy Black Combi



S602/P1658 Patent Black Combi



N302/Y1615 Nubuck Black Combi

Options

Closure Lace or Velcro

Sizes 22 - 45 **Неід**Інт UK & Continental

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles

The following 2 soles are



TR Sole



Spot White



























ARIANA (This style has a medial zip. Rubber sole as standard)





TG1902/L1602 Fabric Black Multicolour



Y1915/L1602 Fantasy Purple



S1997/L1602 Fantasy Blue

Options

Closure Lace or Velcro

Sizes 22 - 45 **Неід**нт UK & Continental

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles

The following 2 soles



TR Sole



Spot White

PU Sole

















Rubber Sole

Black Brown White









L1601 Nappa White

Options

Closure Lace or Velcro

Sizes 22 - 45

Heiqhт UK Only

LASTS & STIFFENERS Standard & Stability 6 widths

Orthosis 3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole





Rubber Welt Sole











N302/L1602 Nubuck Black Combi



V1430/P497 Nappa Blue Combi



R552/AR1920 Waxed Leather Black Combi

Closure Lace or Velcro

Sizes 22 - 45 **Неіqh**т UK Only

Lasts & Stiffeners
Standard & Stability Orthosis 6 widths

3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole






















J846/TR1932 Polished Leather Black Combi



R552/P452 Waxed Leather Black Combi

Closure Lace or Velcro

Sizes 22 - 45 Height UK Only

Standard & Stability 6 widths

Lasts & Stiffeners tability Orthosis 3 widths

Anti-Varus 5 widths

Soles

The following 2 soles



TR Sole





PU Sole

























L1977/AR1977 Nappa Green Combi

M1932/S607 Full Grain Aniline Red Combi



P452/R552 Split Suede Black Combi



V1430/AR1920 Nappa Blue Combi



M1932/NM1902 Full Grain Aniline Grey Combi

Options •

Closure Lace or Velcro

Sizes 22 - 45 **Hεiqh**τ UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

5 widths

Soles =

The following 2 soles



TR Sole





PU Sole























Dakota (PU sole as standard)



Y1955 Fantasy Beige





Y1507 Fantasy Pink



Y1508 Fantasy Black



Y1925 Leopard Fantasy Pink



R1955 Waxed Leather Brown



R552 Waxed Leather Black

Options |

Closure Lace or Velcro

Sizes 22 - 45 **H**eiqhт UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles are



TR Sole



Spot White

PU Sole





Sports Sole









Rubber Sole



















N302/Y1982 Fantasy Black Combi



N328/S1997 Fantasy Blue Combi

Closure Velcro Only

Sizes 22 - 45 Height UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS Orthosis 3 widths

Anti-Varus 5 widths

Soles

The following 2 soles are



TR Sole





PU Sole

















Patrick (TR sole as standard)





R552/N302 Waxed Leather Black Combi



V1430/P497 Nappa Blue Combi



R1955/P407 Waxed Leather Brown Combi

Options

Closure Velcro Only

Sizes 22 - 45 **H**eiqhт UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

5 widths

Soles =

The following 2 soles



TR Sole





PU Sole





Sports Sole









Rubber Sole









(This style has lace closure as standard. Velcro will incur a surcharge)





M1503/P497 Full Grain Aniline Blue Combi



M1909/P1905 Full Grain Aniline Brown Combi



NM1902/P452 Nappa Black Combi

Options

Closure Lace or Velcro

Sizes 22 - 45

Неіднт UK Only

Lasts & Stiffeners Standard & Stability 6 widths

Orthosis 3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole

















L1977/Y1634 Nappa Green Combi



M1903/Y1634 Full Grain Aniline Navy Combi



R1955/Y1634 Waxed Leather Brown Combi

U	PTIONS

Closure Lace or Velcro

Sizes 22 - 45 **Неіqh**т UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

5 widths

Soles =

The following 2 soles



TR Sole





PU Sole





Sports Sole





















R502 Nappa Black *embroidery will incur a surcharge





F1602 Polished Leather Black



L1676 Nappa Blue



TG1902 Fabric Black Multicolour



S606 Patent Fuchsia



S602 Patent Black

\sim		
• ()	DTI	ONS
_	Ь	0113

Closure Lace or Velcro

Sizes 22 - 45 **Н**еі**q**Һт UK & Continental

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles

The following 2 soles are



TR Sole





PU Sole











Sports Sole



Beige

Black Brown White



















- Stops zips from sliding open
- Helps to protect clothes from catching
- Available in the following styles: Kate, Nathan, Sofia, Jessica, Ariana, Tyler, Jean, Jordan & Levi



L1672/S1997 Fantasy Blue Combi



N325/Y1925 Fantasy Pink Combi



L1672/Y1932 Fantasy Black Combi

Options

Closure Lace or Velcro

Sizes 22 - 45

Height UK Only

LASTS & STIFFENERS Standard & Stability 6 widths

<u>Orthosis</u> 3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole











TopGun (This style has lace closure and Sports sole as standard. Velcro will incur a surcharge)





N302/TR1932 Nubuck Black Combi



V1420/P1653 Nappa Navy Combi



R502/P452 Nappa Black Combi

Options

Closure Lace or Velcro

Sizes 22 - 45 **Н**еі**q**Һт UK & Continental

LASTS & STIFFENERS Standard & Stability 6 widths

Orthosis 3 widths

5 widths

Soles

The following 2 soles



TR Sole





PU Sole















Black Brown White

















R507/VS1926 Nappa Fuchsia Combi



R576/L1601 Waxed Leather Blue Combi



N302/TS1976 Nubuck Black Combi

Options

Closure Lace, Velcro or Boa

Sizes 22 - 45 **Неід**Інт UK Only

Standard & Stability 6 widths

Lasts & Stiffeners **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles



TR Sole





PU Sole







Sports Sole





















Goliath (Sports sole as standard)





R502/N302 Nappa Black Combi



V1430/NA1976 Nappa Blue Combi



R577/P496 Waxed Leather Fuchsia Combi

Options •

Closure Velcro Only

Sizes 22 - 45 **H**eiqhт UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

5 widths

Soles =

The following 2 soles



TR Sole





PU Sole































M1401/Y1955 Fantasy Beige Combi



M1401/M1503 Full Grain Aniline Blue Combi



M1401/M1452 Full Grain Aniline Black Combi

Options •

Closure Lace or Velcro

Sizes 22 - 45

НеідІнт UK Only Lasts & Stiffeners Standard & Stability Orthosis 6 widths

3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole















N328/L1601 Nubuck Navy Combi



N1627/L1601 Nubuck Red Combi



N325/L1601 Nubuck Pink Combi

Options

Closure Lace or Velcro

Sizes 22 - 45 **НеіqҺ**т UK Only

LASTS & STIFFENERS Standard & Stability 6 widths

Orthosis 3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole













M1909 Full Grain Aniline Brown





M1903 Full Grain Aniline Navy



Y1932 Fantasy Black



M1452 Full Grain Aniline Black



M1932 Full Grain Aniline Grey



M1927 Full Grain Aniline Red

Options

Closure Lace or Velcro

Sizes 22 - 45 **Н**е**і**q**h**т UK and Continental

LASTS & STIFFENERS Standard & Stability 6 widths

Orthosis 3 widths

Anti-Varus 5 widths

Soles =

The following sole is



TR Welt Sole













Closure System

"Boa Closure" is a unique and stylish closure system available with Nimco footwear!

The Boa system uses steel lace, nylon guides and a mechanical reel that was designed to address deficiencies in traditional closures and offers many advantages when fitted with Nimco footwear.

These include:

- Increased independance for children who have difficulty tying laces allowing easy one handed adjustment with a quick turn of the dial.
- Custom comfort with smooth, even closure and no pressure points.
- Powerful Closure with locked laces staying that way
- Very strong laces made from aircraft grade stainless steel
- Micro adjustability with many degrees of fine tuning
- Extremely light weight and able to shed water, mud and ice!



Please see below styles that are available with the Boa Closure System:



Bambi



Shiny



Snoopy



Cool







Easy to put on and designed specially for use with orthoses







L1672/P402 Nappa Black Combi



N325/Y1925 Nubuck Pink Combi



Y1972/S602 Fantasy Black Combi

Options

Closure **BOA Only**

Sizes 24 - 45 **Hεiqh**τ UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

Anti-Varus 5 widths

Soles =

The following 2 soles



TR Sole





PU Sole





Sports Sole













Rubber Sole







WATER RESISTANT FOOTWEAR

This water resistant style has been designed specifically for people who require all the benefits of Nimco Made4You footwear but also need extra protection for wet weather conditions.

This style is made with special leathers that will not be easily damaged by water. To further increase your protection from rain and snow, you can include NimcoTex Lining* on your order; this membrane is put in between the leather uppers and the lining to further increase the water resistance.

* there is a surchage for NimcoTex Lining





BERG (Rubber sole as standard)









WP590/WP594 Waterproof Leather Brown Combi

Options

Closure Lace or Velcro

Sizes 22 - 45 **Неід**нт UK Only

Standard & Stability 6 widths

LASTS & STIFFENERS **Orthosis** 3 widths

5 widths

Soles =

The following sole is



PU Sole









Rubber Sole



Adaptions

We are pleased to be able to offer a full range of adaptions that can be added to your footwear. As we make our footwear to order, we are able to add many of our adaptions at the time of manufacture. This allows us to maintain a very high quality of product. Below are the adaptions that we currently offer:

- B1 Outside compensating raise up to 25mm with a rand, per shoe
- B2 Each subsequent 25mm or part thereof
- B3 Outside compensating raise up to 25mm without a rand, per shoe
- B4 Each subsequent 25mm or part thereof
- B5 Raising the heel in a single item of footwear up to 25mm
- B6 Lowering the heel
- B7 Addition of heel seat or valgus/varus wedge fitted to shoe, per shoe
- B8 Stirrup, heel, toe raising of caliper retaining strap
- B9 Additional padding in a single item of footwear
- B10 Lower the upper beneath the ankle on one side, per shoe
- B13 Convert to velcro fastening, refacing where necessary, per strap
- B14 Addition of lambswool lining or similar, per shoe
- B15 Addition of high elongated stiffener to the medial/lateral aspect of the upper, per shoe
- B16 Addition of a metatarsal bar fitted to the sole, per shoe
- B17 Inclusion of a rocker sole, per shoe
- B18 Inclusion of a wedge fitted into or onto the sole, per shoe
- B19 Inclusion of a floated sole, per shoe
- B20 Inclusion of a sole strip or plate in composite material, per shoe
- B21 Inclusion of a sole strip or plate in non-composite material, per shoe
- B22 Addition of a rubber or steel tip to the sole or heel, per shoe
- B23 Addition of a tip with a flange to protect the welt, per shoe
- B24 Addition of an external protective toe cap. per shoe
- Addition of a full or half elongated, oblique, floated heel or heel wedge or torque heel, per shoe (Combinations of these items will be charged as multiples of B25)
- B30 Inclusion of a socket, to one or both sides of the heel, without front or backstops, per shoe
- B31 Inclusion of a socket, to one or both sides of the heel, with front or backstops, per shoe
- B32 Inclusion of a toe raising/torsion socket, per shoe
- B33 Addition of a lined and padded T or Y strap, per shoe
- B34 Addition of a lined and padded double decker T or Y strap, per shoe
- Inclusion of a full length or three quarter non moulded inlay in a non-shock absorbent material. bontex, plastazote or leather, covered in a suitable material if required, per shoe
- B40b As B40 but manufactered in shock absorbant material: poron, etc
- B40c As B40 but manufactered in silicone
- Inclusion of a metatarsul or valgus support or heel pad fitted to B40, each
- B42 Inclusion of a metatarsul or valgus support or heel pad fitted directly to the shoe, per shoe
- B43 Inclusion of a heel elevator in cork, HDP or EVA with suitable covering, each
- B50 Inclusion of an inlay with filler to compensate for the amputation of some or all of the toes, per shoe
- B51 Inclusion of an inlay with filler to compensate for the amputation of the forefoot, per shoe
- B54A Functional Foot Orthosis made to cast composite material
- B54C Functional Foot Orthosis or Total Contact inlay made to cast
- B55A Covering for an inlay single density
- B55B Covering for an inlay dual density
- B56 Functional Foot Orthosis posting
- B57 Stretching a single item of footwear

FOOTWEAR REDAIRS:

- C12 Refacing including eyelets, ski hooks per pair
- C13 Replace eyelets, ski hooks or D rings per pair
- C14 Renew laces per pair
- C19 Relining the back of a single item of footwear
- C21 Soling & heeling in synthetic material or rubber, per pair
- C22 Soling & heeling in commando units per pair
- C27 Renew stick on soles per pair



Nimco Stock Inlays

Anti-Pronation Inlay



The anti-pronation foot orthosis aims to resist excessive pronatory moments. It features a medium density EVA base, alcantara top cover and a 3 degrees intrinsic medial heel post.

Neurological Inlay



The neurological foot orthosis incorporates a deep heel cup, peroneal notch, enhanced transverse arch and a 2nd and 5th toe crest with the aim to assist in tone reduction.

LAST Additions

You can specify the following additions to the last when ordering your footwear:

- Extra joint width
- Extra depth (either localised or throughout)
- Extra circumference to the joint or instep
- Changes to the top line circumference
- Hallux Valgus accommodation
- Extra width on base of 5th
- Extra width on the medial side mid foot
- Hammertoes accommodation
- Haglund module

If you would like to add any of these modifications, you can detail them on your order and they will be included when your footwear is manufactured to your specifications.



Pattern Changes

To supplement the range of last additions that we can offer, it is possible to request changes to the pattern of the shoe upper. This allows you to request functional changes such as different opening or cosmetic changes to the design of the shoe to increase compliance.

You can request the following pattern changes to the shoe upper when ordering your footwear.

- Lower opening to make it easier to put shoes on rigid feet or clawing toes
- Addition of a Zip closure to avoid undoing a lace closure each time you need to take your shoes off
- Changes to the design (other than colour) such as extra decorative stripes or rear tabs

If you would like to add any of these modifications, you can detail them on your order and they will be included when your footwear is manufactured to your specifications.



Bespoke To Measures

When your footwear requirements fall completely outside our standard size and width measurements, we offer a fully bespoke to measurement service.

If you would like to order any of our styles based on your custom measures or drafts, please supply the following additional information with your order:

- Size: This can be in mm or expressed as a UK or European shoe size;
- Joint width: This width will ensure that the forefoot and/orthosis will be correctly accommodated;
- Joint circumference: This Circumference allows you to specify the volume required for the forefoot. You can also specify the toe box depth to further tailor your footwear;
- Instep circumference: This measurement ensures that the facings of the shoe will not be too close together or too far apart when closed against the foot;
- Long heel circumference: Providing this measurement will make sure that there is sufficient volume in the rear of the shoe to allow for the rearfoot and ankle or an AFO.

Your order can be processed as a trial fit if necessary to allow further modifications after fitting.



Bespoke To Cast

For more complex foot shapes, we can manufacture any of our styles to a cast of your patient's feet.

In order to create your bespoke to cast order, please send us the patient's foot cast along with the specifications you require the footwear to adhere to. To further improve the accuracy of the cast, you can include a foot outline draft to confirm the correct lengths of the feet.

It is recommended that you have all made to cast orders made to trial fit to ensure that all areas of the feet or AFO are properly accommodated both in sitting and weight bearing positions. Once you have approved the trial, we will apply the sole and finish your footwear to your exact specifications.

Please Note:

Due to the personalised nature of made to measure and made to cast orders, they are treated as non returnable.



GENERAL TERMS

Ordering

To place an order please post, email or fax us with your requirements. If you prefer to use product codes to describe your footwear, please contact us and we will be happy to send you our latest order code list.

Delivery

Custom Made Footwear is made to order and so we are able to provide very accurate delivery times. Once we have received your order it will be processed and sent electronically to Nimco for manufacture. On receipt of your order we will deliver to you in 12 - 14 working days. This delivery timescale is very reliable so you are able to make your next patient appointment with confidence. If a delay does occur with your order, we will always contact you and explain the delay.

Split Pairs

Split pairs are processed in the same way as the custom made range unless agreed otherwise. Once your order is received it will attract the same reliable delivery scale as the custom made range. All split pair orders attract a 50% surcharge and can be split in any way (last. size. stiffener).

RETURNS

When you receive your order we want you to be happy with every aspect of your footwear. However if you need to return the footwear for any reason. We will be happy to receive them under the following conditions:

- Your footwear must be unworn and be less than 30 days from date of delivery.
- Most custom made boots are returnable, however they attract a 25% restocking fee.
- Unfortunately, we cannot accept returns of sandal styles, i.e. Goliath and Surf.
- Footwear which has been ordered as a split pair are also non-returnable.
- Finally, any order which has been adapted in any way or utilises non catalogue colours or styles can not be returned.

Your statutory rights are not affected.

Payment And Delivery

Payment terms and delivery tariff information are available on request as this is dependent on your geographical location and pre agreed credit terms. If you have any other queries. please don't hesitate to contact us and we will be more than happy to assist you.













Standard Lasts Measurements

This table illustrates the measures for the standard lasts.

The Length⁽¹⁾, Joint Width⁽²⁾ and Instep circumference⁽³⁾ positions are illustrated below for reference. All sizes are with inlay inserted into shoe.

Shoe* TN1 -X		Narrow	w TN2 -Narrow		TN4 -Medium		TN6 - Wide		TN8 - X Wide		TN10 - XX Wide		
Size	Length(1)	Joint(2)	Instep(3)	Joint(2)	Instep(3)	Joint(2)	Instep(3)	Joint(2)	Instep(3)	Joint ⁽²⁾	Instep(3)	Joint(2)	Instep(3)
16	105	43	116	45	121	49	131	53	140	56	150	60	160
17	112	44	120	46	125	50	135	54	144	57	154	61	164
18	118	46	124	48	129	52	138	56	148	59	158	63	168
19	125	48	128	50	133	54	143	58	152	61	162	65	172
20	132	50	132	52	137	56	147	60	156	63	166	67	176
21	138	52	136	54	141	58	151	62	160	65	170	69	180
22	145	53	140	55	145	59	155	63	164	66	174	70	184
23	152	54	144	56	149	60	159	64	168	67	178	71	188
24	158	56	149	58	153	62	163	66	172	69	182	73	192
25	165	57	152	59	157	63	167	67	176	70	186	74	196
26	172	59	156	61	161	65	171	69	180	72	190	76	200
27	178	60	160	62	165	66	175	70	184	73	194	77	204
28	185	61	164	63	169	67	179	71	188	74	198	78	208
29	192	63	168	65	173	69	183	73	192	76	202	80	212
30	198	64	171	66	176	70	186	74	195	77	206	81	216
31	205	65	175	67	180	71	190	75	199	78	210	82	220
32	212	67	178	69	183	73	193	77	202	80	214	84	224
33	218	69	182	71	187	75	197	79	206	82	218	86	228
34	225	70	186	72	191	76	201	80	210	83	222	87	232
35	231	71	190	73	195	77	205	81	214	84	226	88	236
36	238	73	194	75	199	79	209	83	218	86	230	90	240
37	245	74	198	76	203	80	213	84	222	87	234	91	244
38	251	76	202	78	207	82	217	86	226	89	238	93	248
39	258	77	206	79	211	83	221	87	230	90	242	94	252
40	265	79	211	81	216	85	226	89	235	92	246	96	256
41	271	81	215	83	220	87	230	91	239	94	250	98	260
42	278	82	219	84	224	88	234	92	243	95	254	99	264
43	285	84	223	86	228	90	238	94	247	97	258	101	268
44	291	86	227	88	232	92	242	95	251	98	262	102	272
45	298	88	231	90	236	94	246	97	255	100	266	104	276

All dimensions are in millimeters

The AV Lasts TN2/AV, TN4/AV, TN6/AV, TN8/AV and TN10/AV are comparable to the standard last measurements but with a different forefoot alignment



Measurements:

1 - Length 2 - Joint Width

3 - Instep Circumference





^{*}Please ensure that desired toe room is added to the length before choosing shoe size

Orthoses Lasts Measurements

This table illustrates the measures for the orthoses lasts. The TN8 and TN10 lasts are the same as the standard lasts but with reduced toe spring. The TN10xls last is the same as the TN10 but with extra volume by increasing the long heel circumference.

The Length(1), Joint Width(2), Instep circumference(3) and Long Heel circumference(4) positions are illustrated below for reference. All sizes are with inlay inserted into shoe.

Shoe*		TN	8/ort - X 1	Wide	TN10	ort - XX	Wide	TN10xls/ort - Extra Volume			
Size	Length(1)	Joint(2)	Instep ⁽³⁾	Long Heel ⁽⁴⁾	Joint ⁽²⁾	Instep(3)	Long Heel ⁽⁴⁾	Joint ⁽²⁾	Instep(3)	Long Heel ⁽⁴⁾	
16	105	56	150	175	60	160	180	60	160	189	
17	112	57	154	183	61	164	188	61	164	197	
18	118	59	158	192	63	168	197	63	168	206	
19	125	61	162	200	65	172	205	65	172	214	
20	132	63	166	208	67	176	213	67	176	222	
21	138	65	170	217	69	180	222	69	180	231	
22	145	66	174	225	70	184	230	70	184	239	
23	152	67	178	234	71	188	239	71	188	248	
24	158	69	182	242	73	192	247	73	192	256	
25	165	70	186	251	74	196	256	74	196	265	
26	172	72	190	259	76	200	264	76	200	273	
27	178	73	194	268	77	204	273	77	204	282	
28	185	74	198	276	78	208	281	78	208	290	
29	192	76	202	285	80	212	290	80	212	299	
30	198	77	206	294	81	216	299	81	216	308	
31	205	78	210	302	82	220	307	82	220	316	
32	212	80	214	310	84	224	315	84	224	324	
33	218	82	218	319	86	228	324	86	228	333	
34	225	83	222	328	87	232	333	87	232	342	
35	231	84	226	336	88	236	341	88	236	350	
36	238	86	230	345	90	240	350	90	240	359	
37	245	87	234	353	91	244	358	91	244	367	
38	251	89	238	362	93	248	367	93	248	376	
39	258	90	242	370	94	252	375	94	252	384	
40	265	92	246	378	96	256	383	96	256	392	
41	271	94	250	387	98	260	392	98	260	401	
42	278	95	254	396	99	264	401	99	264	410	
43	285	97	258	404	101	268	409	101	268	418	
44	291	98	262	412	102	272	417	102	272	426	
45	298	100	266	421	104	276	426	104	276	435	

All dimensions are in millimeters

*Please ensure that desired toe room is added to the length before choosing shoe size

The AV Lasts TN2/AV, TN4/AV, TN6/AV, TN8/AV and TN10/AV are comparable to the standard last measurements but with a different forefoot alignment



Measurements:

1 - Length

2 - Joint Width

3 - Instep Circumference







DENOVO HEALTHCARE LTd Unit 6a Cheddar Business Park

Unit 6a Cheddar Business Park Wedmore Road Cheddar BS27 3EB

Tel: 01934 808 416 Fax: 01934 808 405

E-mail: sales@denovohealthcare.com Website: www.denovohealthcare.com

