

Adult Diabetic Footwear Range

2019/2020 Collection

Now Available
-NEW-
in 2 widths



est. 1904
NIMCO
MADE 4 YOU



Contents

Collection Characteristics & Specifications	2
Womens	4
Mens	12
Insoles Characteristics & Specifications	20
General Terms & Conditions	21
Measurement Charts	22

Brand New Diabetic Footwear Collection

In this collection you will find 13 different styles. All of the shoes are fashionable and comfortable. 6 styles have our new sporty PU sole which matches the current trend in the shoe industry; 3 styles are made with a light weight EVA sole which makes them a great recommendation if you are looking to use them as indoor shoes.

Characteristics of this new collection:

- A wide and diverse range of models - all models are stock items.
- All models are produced using our approved ladies last K/M and men's last K/M: nature shaped front part providing sufficient space for toes – narrow heel part providing a perfect fitting which avoids heel slippage.
- All models are produced with “Xsensible Inside Stretch Leather” – this provides greater flexibility to accommodate for changes in foot volume, eg oedema or prominent bunions.
- Special diabetic lining material designed for diabetic patients - soft, comfortable and durable which is also breathable, washable and skin friendly
- Shoes include 2 inlays: a 5mm cork insole and a 5mm cork-leather insole - both are removable and the maximum depth for individual insoles is 10mm.
- Additional reinforced insoles are available upon request, not permanently integrated in the shoe.
- All shoes are available in UK ladies sizes 3-9 and men sizes 6-12. (including ½ sizes)
- Outsoles: PU outsoles and light weight EVA. Both easy to adapt.





Xsensible Stretch Leather

The advantages and characteristics:

- This feature ensures that the leather stretches 70% in width and also 30% in length.
- The top layer is made of genuine leather and the bottom layers are made of Lycra constructed with a natural foam material in between.
- All footwear in this collection is made with Xsensible Stretch Leather technology which is extremely stretchable and adapts individually to your foot.
- Xsensible Stretch leather has a unique property of being breathable and providing an optimal foot temperature

Indications

Diabetic polyneuropathy with pronounced Neuropathy (loss of sensitivity) and more pronounced Angiopathy (circulatory disorders), associated with already existing or healed ulcers and/or other mechanical irritations that lead to regional pressure peaks.



DeNovo Healthcare



Quality • Service • Choice

Find us:



@denovohealthcareltd_uk



@DeNovoHealthcareLtd



@denovoUK



W
O
M
E
N
S





Sand nubuck/leather with sand stretch leather



Black leather/nubuck/patent with black stretch nubuck

Sizes Available: 3 to 9 (Also available in 1/2 sizes)

Last: K & M

Sole: PU Sport

Matching sole stiffeners: see page 18



Grey fantasy suede with grey stretch leather

Sizes Available: 3 to 9 (Also available in 1/2 sizes)

Last: K & M

Sole: PU Sport

Matching sole stiffeners: see page 18



Black leather with black stretch leather
PU Sole



Anthracite leather with anthracite stretch leather
EVA Sole

Sizes Available: 3 to 9 (Also available in 1/2 sizes)

Last: K & M

Sole: PU Sport / Lightweight EVA

Matching sole stiffeners: see page 18





Taupe fantasy leather with taupe stretch leather



Black leather with black stretch leather

Sizes Available: 3 to 9 (Also available in 1/2 sizes)

Last: K & M

Sole: PU Sport

Matching sole stiffeners: see page 18



Black nubuck with black stretch leather



Style constructed with medial zip for easy opening

Sizes Available: 3 to 9 (Also available in 1/2 sizes)

Last: K & M

Sole: PU Sport


Matching sole stiffeners: see page 18



Pashmina Croco leather with taupe stretch leather



Navy metal leather/nubuck with navy stretch leather

 Style constructed with lateral zip for easy opening

Sizes Available: 3 to 9 (Also available in 1/2 sizes)

Last: K & M

Sole: PU Sport

Matching sole stiffeners: see page 18

M
E
N
S



Brave Lace



Grey nubuck with grey stretch leather

Sizes Available: 6 to 12 (Also available in 1/2 sizes)

Last: K & M

Sole: Sporty lightweight EVA Sole

Matching sole stiffeners: see page 18



Black nubuck with black stretch leather

Sizes Available: 6 to 12 (Also available in 1/2 sizes)

Last: K & M

Sole: Sporty lightweight EVA Sole

Matching sole stiffeners: see page 18



Dark blue nubuck & black leather with navy stretch leather

Sizes Available: 6 to 12 (Also available in 1/2 sizes)

Last: K & M

Sole: Sporty lightweight EVA Sole

Matching sole stiffeners: see page 18

Sailor Lace



Grey leather with grey stretch leather

Sizes Available: 6 to 12 (Also available in 1/2 sizes)

Last: K & M

Sole: PU Sole

Matching sole stiffeners: see page 18

Sailor Velcro



Brown leather with brown stretch leather



Sizes Available: 6 to 12 (Also available in 1/2 sizes)

Last: K & M

Sole: PU Sole

Matching sole stiffeners: see page 18



Black leather with black stretch leather

Sizes Available: 6 to 12 (Also available in 1/2 sizes)

Last: K & M

Sole: PU Sole

Matching sole stiffeners: see page 18



Black leather with black stretch leather

Sizes Available: 6 to 12 (Also available in 1/2 sizes)


Last: K & M

Sole: PU Sole

Matching sole stiffeners: see page 18

Insoles Characteristics & Specifications

All of our diabetic stock footwear include two removable separate insoles:



By removing both insoles, you gain an extra 10mm depth and 20mm in circumference.

Having removable insoles also makes it possible for a bespoke orthosis to be worn with these shoes comfortably.



Sole reinforcements **VIPREX®**

As an addition it is also possible to order a perfectly fitting sole stiffener made of innovative glass fibre polypropylene. This latest development is extremely thin but has the highest rigidity within the shoe. It is a simple and health-friendly component. Thickness: 1mm

Available sizes:

- | | |
|------------|---------------------------------------|
| Ladies I | for shoe size 3 - 4.5 (code 9114/1) |
| Ladies II | for shoe size 5 - 7.5 (code 9114/2) |
| Ladies III | for shoe size 8 - 9.5 (code 9114/3) |
| Men's I | for shoe size 6 - 7.5 (code 9115/1) |
| Men's II | for shoe size 8 - 10.5 (code 9115/2) |
| Men III | for shoe size 11 - 12.5 (code 9115/3) |



Tip: To achieve an even firmer sole stiffness on the ladies shoes use a mens stiffener as it has a higher rigidity.

General Terms

Ordering

To place an order, please post, email or fax us with your requirements. If you prefer to use product codes to describe your footwear, please contact us and we will be happy to send you our latest order code list.

Delivery

Our Diabetic Stock Footwear is part of a stock system. Once we have received your order it will be processed and dispatched within 2 to 3 working days (subject to stock). If a delay does occur with your order, we will always contact you, explain the delay and provide a new delivery date.

Returns

When you receive your order we want you to be happy with every aspect of your footwear. However if you need to return the footwear for any reason. We will be happy to receive them under the following conditions:

- Your footwear must be unworn and be less than 30 days from date of delivery.
- Finally, any order which has been adapted in any way cannot be returned unless faulty.

Your statutory rights are not affected.

Payment And Delivery

Payment terms and delivery tariff information are available on request as this is dependent on your geographical location and pre agreed credit terms. If you have any other queries. please don't hesitate to contact us and we will be more than happy to assist you.

VAT

VAT will be applied at the ruling rate. Where VAT exemption is applicable, an Exemption Certificate or accepted documentation must be supplied at the point of ordering.

PAYMENT TERMS

Payment terms are strictly 30 days from date of invoice for account holders. BACS, Cheque, MasterCard or Visa may also be used to make payment. Please allow 7 working days for cheque clearance before goods are dispatched. Please make all cheques payable to DeNovo Healthcare Ltd. We reserve the right to charge interest on overdue credit accounts and severely overdue accounts may be placed on stop. We will always make every attempt to contact you to discuss your account before this action occurs.

Please refer to our Pricelist for full Terms & Conditions, these can be requested via phone, email or letter.

Ladies Measurement Charts

Last K (normal)

Size	1	2	3	4	5	6	7	8
3	235	88	224	228	304	20	23	15
3.5	239	89	228	231	309	21	24	15
4	243	90	231	234	314	21	24	16
4.5	247	91	234	237	320	21	24	16
5	252	92	237	241	325	22	25	16
5.5	256	93	240	244	330	22	25	16
6	260	94	243	247	335	22	25	16
6.5	264	95	247	250	340	22	25	17
7	268	96	250	253	345	23	26	17
7.5	273	97	253	256	350	23	26	17
8	277	98	256	260	355	23	26	17
8.5	281	99	259	263	360	24	27	17
9	285	100	262	266	365	24	27	18

Last M (wide)

Size	1	2	3	4	5	6	7	8
3	235	91	234	238	314	21	23	15
3.5	239	92	238	241	319	22	24	15
4	243	93	241	244	324	22	24	16
4.5	247	94	244	247	330	22	24	16
5	252	95	247	251	335	23	25	16
5.5	256	96	250	254	340	23	25	16
6	260	97	253	257	345	23	25	16
6.5	264	98	257	260	350	23	25	17
7	268	99	260	263	355	24	26	17
7.5	273	101	263	266	360	24	26	17
8	277	102	266	270	365	24	26	17
8.5	281	103	269	273	370	25	27	17
9	285	104	272	276	375	25	27	18

All of the measurements are in millimeters. These measurements are with the inlay in the shoes, you can achieve up to 10mm extra depth and/or 20mm extra circumference by removing both insoles.

Mens Measurement Charts



Last K (normal)

Size	1	2	3	4	5	6	7	8
6	262	99	240	246	344	23	23	9
6.5	266	100	244	250	349	24	24	9
7	271	101	247	253	354	24	24	10
7.5	275	102	250	256	360	24	24	10
8	279	103	253	259	365	25	25	10
8.5	283	104	256	262	370	25	25	10
9	288	105	259	265	375	25	25	10
9.5	292	106	263	269	380	25	25	10
10	296	107	266	272	385	26	26	10
10.5	300	108	269	275	390	26	26	10
11	304	109	272	278	395	26	26	11
11.5	309	110	275	281	400	27	27	11
12	313	112	278	285	405	27	27	11

Last M (wide)

Size	1	2	3	4	5	6	7	8
6	262	102	250	256	354	24	23	9
6.5	266	103	254	260	359	25	24	9
7	271	104	257	263	364	25	24	10
7.5	275	105	260	266	370	25	24	10
8	279	106	263	269	375	26	25	10
8.5	283	107	266	272	380	26	25	10
9	288	109	269	275	385	26	25	10
9.5	292	110	273	279	390	26	25	10
10	296	111	276	282	395	27	26	10
10.5	300	112	279	285	400	27	26	10
11	304	113	282	288	405	27	26	11
11.5	309	114	285	291	410	28	27	11
12	313	115	288	294	415	28	27	11

- 1 Foot Length
- 2 Joint Width
- 3 Joint Circumference
- 4 Instep Circumference
- 5 Long Heel
- 6 Toe Height
- 7 Heel Height
- 8 Toe Spring





DeNovo Healthcare Ltd

6a Cheddar Business Park
Wedmore Road
Cheddar, Somerset BS27 3EB

Tel: 01934 808 416

Fax: 01934 808 405

E-mail: sales@denovohealthcare.com

Website: www.denovohealthcare.com

Proud suppliers for:

